

OrthoSave® Classic 2 Compression Therapies *Patient User Guide*

Basic Purpose of the OrthoSave® Classic 2 Compression Therapy

To provide controlled compression therapy for the reduction of swelling and pain, as well as increasing joint range of motion. In addition, the OrthoSave® system provides mechanical blood clot prevention.

Your Daily Prescription (default prescription in parentheses)

Use for _____(30) minutes with 1 hour off between sessions _____(3) times per day for _____(16) days.

Basic Patient Guide for Use of the OrthoSave® Classic 2

- Step 1. Place all wraps on appropriate body parts.
- Step 2. Elevate affected body parts if possible or reasonable.
- Step 3. Turn machine on at flip switch in the back.
- Step 4. Press DVT settings button for foot or calf if your rep demonstrated the need to do so.
- Step 5. When session is completed, turn off at flip switch in the back.

Basic Troubleshooting

- ***If needing to stop in mid therapy:***
Flip “OFF” switch. Therapy will stop.
Disconnect wraps from your body.
- ***If flow of therapy blocked:***
Check the lines adjoining the wrap to the machine for proper flow. If the therapy is not flowing properly, make sure the lines are not kinked then vent the system while the therapy is running to enable trapped air pressure to escape by opening the top reservoir. Shake the wraps to loosen air pockets in the lines. Screw cap back on and resume normal therapy. If you are using the therapy for a foot, reposition your foot to make sure there is space between the wrap and where you are resting it for air to overcome the pressure of your leg against the furniture’s surface. Restart therapy.



For all other questions: Contact _____ at _____.
(Rep Name)